

Silicone Scar-care Dressing

Introduction

Far Eastern Group used medical-grade silicone gel to develop self-adhesive scar-care dressing (Fig. 1). The dressing has the abilities of moisturizing and pressing which can maintain the scar under the pressure and in the moisture environment, and then effectively inhibits scar formation to improve the appearance of scars. The dressing with nylon layer in skin color can conceal the scar and block out 98% UV rays to avoid the pigmentation due to sunlight.

The clinical evidence shows that the silicone scar-care dressing can effectively prevent the scar formation using on close wounds (Fig. 2). There are obvious improvements on the wounds due to surgery, accident, trauma, keloids and hypertrophic scars.




Fig. 1 Silicone Scar-care Dressing

Clinical Result:

Female/Surgical Wound/Period: 2008.11.25-2009.06.01



Fig. 2 Pictures of clinical test result

Properties:

- ✓ self-adhesiveness
- ✓ no residue, no degumming
- ✓ flexible and perfectly fit on your scar
- ✓ pass SGS UV block test (skin color)
- ✓ washable by neutral detergent and can be reused for 32 times
- ✓ medical-grade silicone ingredients
- ✓ high biocompatibility

Applications

- scars caused by surgery (e.g., caesarean section, minimal invasive surgery)
- scars caused by trauma (e.g., puncture, knife wound)
- scars caused by burns
- keloids and hypertrophic scars