


Hydrocolloid dressings

Introduction

Hydrocolloid dressing consists of a hydrocolloid adhesive with an outer transparent covering film impermeable to fluids and bacteria. The inner layer of hydrocolloid adhesive can absorb excess wound exudate, maintain moisture balance in the wound, and ensure an optimal moist wound environment. Hydrocolloid dressing has suitable adhesive property and can effectively protect the wound; furthermore, it resembles the color of human skin, as a result, it can keep a natural appearance.



Product Features

- has suitable adhesive property and can effectively protect the wound
- absorbs excess wound exudate
- provides a self-debridement environment
- occlusive, and can induce angiogenesis by providing an anaerobic environment for the wound.
- waterproof, bacteria proof and moisture permeable
- its color is close to human skin color, as a result, it can keep a natural appearance

Indications

- Superficial wounds
- Partial- and full-thickness dermal ulcers
- Post-operative wounds
- Donor sites