

# Liquid Bandage

## Introduction

**Liquid Bandage** is an alcohol-free liquid dressing which consists of silicone materials. It is applied to the damaged or fragile skin and forms a protective barrier which is clear, flexible, breathable, and waterproof and acts as a protective covering allowing the wound to heal. It remains on skin for several days, and naturally wears off as wound heals.



**Fig. 1 Liquid Bandage**

## Properties

- ✓ Quick-drying liquid forms a waterproof, breathable and hypoallergenic protection film on the skin surface
- ✓ Protect wounds and periwound skin from second damage
- ✓ No irritant ingredients to reduce the painful feeling of patients

- ✓ The protection film naturally wears off with the skin metabolism within 72hrs
- ✓ Pass the SGS biocompatibility tests

## Applications

Liquid Bandage is used in:

- Surgical wound (after removing surgical suture )
- Protecting fragile skin and injury wound
- Lacerations
- Abrasion